



TheraCem®

Dual-Cured, Self-Adhesive Calcium & Fluoride Releasing Cement

Duo-Link Universal™

Dual-Cured Resin Luting Cement

Choice™ 2

Light-Cured Veneer Cement

BisCem®

Dual-Cured, Self-Adhesive Resin Cement

eCement® LC

Light-Cured Veneer Cement for Lithium Disilicate

eCement® DC

Dual-Cured Resin Luting Cement for Lithium Disilicate

BISCO Cements

Curing Mode	Dual Cured	Dual Cured	Light Cured	Dual Cured	Light Cured	Dual Cured
Requires Adhesive		●	●		●	●
Substrates	zirconia, glass ceramics, lithium disilicate, metal, alumina, composite	zirconia, glass ceramics, lithium disilicate, metal, alumina, composite	zirconia, glass ceramics, lithium disilicate, composite	zirconia, glass ceramics, lithium disilicate, metal, alumina, composite	lithium disilicate	lithium disilicate
Inlays/Onlays	●	●		●		●
Crowns/Bridges	●	●		●		●
Endodontic Posts	●	●		●		●
Maryland Bridges		●				●
Veneers			●		●	
Primer Required		●	●	●	●	●
Silane recommended for Glass Ceramics/Porcelain or Lithium Disilicate/Silicate	●	●	●	●	●	●
Best Uses	A self-adhesive cement, ideal for retentive preparations when "bonding" the restoration is NOT indicated. Releases both calcium and fluoride. Contains MDP.	Ideal for less retentive preparations when "bonding" the restoration IS indicated. Ideal for most restorations including CAD/CAM restorations.	Ideal cement for veneer restorations. Choice 2 is specifically formulated for high esthetics with color stability.	Translucent, dual-cured, self-adhesive resin cement. Indicated when an esthetic self-adhesive cement is preferred.	Designed for lithium disilicate veneer restorations (including CAD/CAM veneers).	Designed for lithium disilicate restorations (including CAD/CAM restorations).

INDICATIONS



Consult Instructions For Use

Rx Only

MC-51193RC

1-800-247-3368 • www.bisco.com

