


BISCO

TheraCal PT[®]

Dual-Cured


Dual-Cured Resin-Modified Calcium Silicate Pulpotomy Treatment

1




Perform pulpotomy treatment.

2




Achieve hemostasis.

3



Place TheraCal PT directly in the pulp chamber, and ensure good adaptation to the cavity walls and margins.

4



Light cure TheraCal PT for 10 seconds. Place desired adhesive, base, and/or restoration.



Rx Only

 BISCO, Inc. • 1100 W. Irving Park Road • Schaumburg, IL 60193 • U.S.A. • 1-847-534-6000 • www.bisco.com

 BISICO France • 208, allée de la Coudoulette • 13680 Lançon de Provence • France • 33-4-90-42-92-92



MC-4014TP